

AccuCup[®] Dosage Cups and Spoons



Because Dosing Accuracy is Important

Everyday, patients and consumers rely on dosage cups and spoons to measure out and administer a wide-variety of prescription and over-the-counter liquid medications. But often the clear embossed lines on them prove challenging to read, especially in lowlight conditions, introducing the potential for dosing error.

In response, we created the AccuCup[®] line of dosing solutions with boldly printed, dishwasher-safe, easy to read lines. During production our in-line vision systems inspect each cup and spoon to assure the measurement lines are accurately printed and precisely placed, so you can trust you are getting the correct dose, every time.



**73 YEARS OF PACKAGING
LEADERSHIP—AND THE REST IS HISTORY.**

Visit comar.com or dial 1-866-655-7165 for more info.

AccuCup®

Printed Dosage Cup, Round



- Printed measurement markings for improved visibility
- Markings are in milliliters (mL), in response to regulatory guidance
- Standard blank cup design accommodates custom printing for specific measurements, brand, and dosing instructions
- 360° printing in up to 2 colors
- FDA-compliant, food-grade inks are UV-cured making our cups washable and reusable
- Durable printing on the outside - no contact with medication
- Can be custom-cored for a perfect bottle cap fit

Item Family	Description	Cup Size	Resin	Cup Type
22-0481	Round AccuCup®, Natural	20 mL	PP	Printed
22-1211	Round AccuCup®, Natural	30 mL	PP	Printed
22-0717	Round AccuCup®, Natural	60 mL	PP	Printed

Printed Dosage Spoon



- Printed measurement markings for improved visibility
- Markings are in milliliters (mL), in response to regulatory guidance
- Durable printing on the outside - no contact with medication
- FDA-compliant, food-grade inks are UV-cured making our cups washable and reusable
- Custom printing options available

Item Family	Description	Cup Size	Resin	Cup Type
22-1549	Dosage Spoon, Natural	10 mL	PP	Printed

AccuCup®

Printed Dosage Cup, Square



- Printed measurement markings for improved visibility
- Markings are in milliliters (mL), in response to regulatory guidance
- Standard 20 ml printed cup but in a square form
- Custom imprinting available on two opposing sides in two colors
- FDA-compliant, food-grade inks are UV-cured making our cups washable and reusable
- Durable printing on the outside - no contact with medication
- Fits several 24 mm closures

Item Family	Description	Cup Size	Resin	Cup Type
22-1525	Square AccuCup®, Natural	20 mL	PP	Printed
22-1525	Square AccuCup®, Blue Tinted	20 mL	PP	Printed

Embossed Dosage Cup



- Embossed markings for 5-30 mL dosages
- Markings are in milliliters (mL), in response to regulatory guidance

Item Family	Description	Cup Size	Resin	Cup Type
22-0263	Dosage Cup, Natural	30 mL	PP	Embossed

0280 DOSAGE 04.22

**73 YEARS OF PACKAGING LEADERSHIP—
AND THE REST IS HISTORY.**

Visit comar.com or dial 1-866-655-7165 for more info.